

Table1: KML 12-03

Ni Assay (AR-ICP)	Ni Assay (TD-ICP)	Core Length*	From (m)	To (m)
0.17% Ni	0.20% Ni	250.4m	14.0	264.4
including:				
0.23% Ni	0.27% Ni	122.4m	142.0	264.4

(*Note: all measurements are down-hole core lengths, true width has not yet been calculated)

Table2: KML 12-04

Ni Assay (AR-ICP)	Ni Assay (TD-ICP)	Core Length*	From (m)	To (m)
0.22% Ni	0.25% Ni	114.2m	314.0	428.2
including:				
0.24% Ni	0.27% Ni	72.2m	356.0	428.2
Including				
0.26% Ni	0.29% Ni	11.2m	417.0	428.2

(*Note: all measurements are down-hole core lengths, true width has not yet been calculated)